




**ROCKMED<sup>®</sup>**  
*Industries*

Extraction / Extraktion / Extracción / Estrazione



**Dental Instruments catalogue**



*Precision Craftsmanship*

# PDL Luxating Elevators

Precision and durability combined



## PDL Accu Luxating Set



MDE.PDL.07A Set of 7



MDE.PDL.07B Set of 7



MDE.PDL.01 Distal, mm 4, Micro edge    MDE.PDL.02 Mesial, mm 4, Micro edge    MDE.PDL.03 Distal, mm 2.5, Micro edge    MDE.PDL.04 Mesial, mm 2.5, Micro edge    MDE.PDL.05 mm 2,5 Micro edge, curved    MDE.PDL.06 mm 2,5 Micro edge, straight    MDE.PDL.07 mm 3,4 Straight

## PDL Accu Bernard Spade Set



MDE.PDL.14A Set of 7



MDE.PDL.14B Set of 7

ROCKMED Titanium Blue Luxating Elevator Serrated




Antibacterial titanium coating!  
Easy to clean! Safe!



MDE.PDL.08 mm 3,4, Left    MDE.PDL.09 mm 3,4, Right    MDE.PDL.10mm 3,4, External    MDE.PDL.11 mm 3,4, Internal    MDE.PDL.12 mm 4,5, Internal    MDE.PDL.13 mm 4,5, External    MDE.PDL.14 mm 3,4, Straight



MDS.300.13N Set of 7 PDL Elevator - Titanium Black



MDS.300.13 Set of 7 Individually Packed ( Without Cassette )

## ROCKMED PROXI PERIROTOMES SET OF 7

Set of 7 ROCKMED Proxi Perirotomes, including : Small Straight , Small Curved , Sharpened Spade , Small Mesial and Distal , Large Mesial and Distal.

MDS.300.17X Set of 10 in Cassette Hold 10 7 PDL, 3 Splitting Chisels

MDS.300.14 Set of 7 ( Cassette 10 Hold Single Deck )



## Spade Root Forceps

## Spade Root Elevators 4mm

For clinicians who are interested in the entire, Medflair has created a kit with all instruments for atraumatic extractions.



Set of 8 ( 2 Spade Root Forceps & 6 Spade Root Elevators ) MDE.95.10

Set of 6 ( 6 Spade Root Elevators ) MDE.95.09

MDE.95.01  
UPPER ROOT,  
SERRATED

MDE.95.02  
LOWER ROOT,  
SERRATED

MDE.95.03  
MESIAL

MDE.95.04  
RIGHT

MDE.95.05  
STRAIGHT

MDE.95.06  
CURVED

MDE.95.07  
LEFT

MDE.95.08  
DISTAL

## PERIOTOMES

- used at the beginning of any extraction, flexible to follow the curvature of the root
- feature a sharpened edge along the entire radius, thin enough to apply along the fragile buccal wall
- most conservative ligament incising instruments available, the original, ROCKMED instruments for atraumatic exodontia
- sever the periodontal ligament without expansion to the alveolus
- available in ROCKMED three configurations; Anterior, Posterior and Perpendicular
- inserted vertically into the PDL space, parallel to the long axis of the root
- titanium nitride coated for durability and edge retention; they never need sharpening
- reduce the need for traditional elevation and reduce rotational luxation with extraction forceps

Antibacterial titanium coating !  
Easy to clean ! Safe !



SET OF TWO  
MDE.96.010.011  
PERIOTOME #1  
PERIOTOME #2

PERIOTOME #1 MDE.96.010

TiNi COATED PERIOTOME #1, ANTERIOR,  
STRAIGHT BLADE / SOCKET EXPANDER



PERIOTOME #2 MDE.96.011

TiNi COATED PERIOTOME #2, POSTERIOR,  
45° CONTRA-ANGLED BLADES



PERIOTOME #3 MDE.96.012


TiNi COATED PERIOTOME # 3, PERPENDICULAR BLADES





### HEIDBRINK root apex elevators

Fine, pointed elevators with extremely long working ends for removing fractured root apex.



### MDE.RE.HF9 Set of 3

#### Luxating Elevator/Periotome DUAL EDGE

For Reduced Trauma Extractions with 3mm, 4mm, 5mm

Luxating/Periotome Combination Elevator

Delicate, Sharp Instruments Cut Periodontal Ligament While Preserving Bone Anatomy

Thin, Curved, Tapered Tips Ensure Smooth Widening Of The Alveole

Fine Blade Inserts Between Bone And Root With Minimal Damage To Bone Structure

Alternative To Standard Elevators For Less Traumatic Extractions



### Fingerform LUXATORS

#### MDE.RE.HF Set of 8

-For Reduced Trauma Extractions

-Delicate, Sharp Instruments Cut Periodontal Ligament While Preserving Bone Anatomy


-Thin, Curved, Tapered Tips Ensure Smooth Widening Of The Alveole

-Fine Blade Inserts Between Bone And Root With Minimal Damage To Bone Structure

-Alternative To Standard Elevators For Less Traumatic Extractions

-Set Includes Straight And Curved Tips For Access To Both Anterior And Posterior Teeth

-The sharp and flat working tip for loosening of circular ligament. Subsequently the extraction is possible without hurt of mucous membrane




### Heidbrink Root Tip Elevators





# Luxators / Luxateurs Luxadores / Luxationsinstrumente



Luxators  
Set of 10 Pieces  
**MDE.LUXS**




Luxators  
Set of 10 Pieces  
Micro Serrated  
**MDE.LUXSR**

## ROCKMED LUXATORS



**Set Code # MF-LUX-SET**  
Set of 10 Pieces



Designed to be easily inserted along with root surface to access and cut the periodontal ligament prior to extraction.  
Using luxator will minimize tissue trauma and promote quicker healing.  
A real alternative to traumatic forceps extractions.  
Long slender profile, standard - total length: 17,5 cm  
Original Luxators with Color Coded Identification

Available in Original Polished, Titanium Black, Titanium Blue, Titanium Gold TiNi  
Set of 10 (L-1S, L-2S, L-3S, L-3C, L-3CA, L-3IC, L-5C, L-5S, DE-3, DE-5)



Add - "G" For Gold

# Root Elevators | Root-Aufzüge Ascenseurs de racine | Elevadores de raíz



The Angled Luxators are sharp enough to cut into the PDL, thin enough at the tip so that it can fit under the tooth below the gumline and strong enough to elevate it without breaking. They minimize the time and effort it takes to switch between things like periostomes, elevators, and luxators while also accomplishing an atraumatic extraction, 3mm Right angle, 3mm Left angle, and 5mm Straight tipped.



## KOPP ROOT ELEVATORS



1 KOPP

2 KOPP

3 KOPP

4 KOPP

5 KOPP

6 KOPP



3,6 mm  
Upper roots

3,6 mm  
Upper roots

4,2 mm  
Upper roots

4,2 mm  
Upper roots

4,2 mm  
Lower roots

4,2 mm  
Lower roots

## For Children



#301V Modified Apexo 301  
2.0 mm blade width.  
MDE.RE.54

Apexo #301 Straight  
2.5 mm blade width.  
MDE.RE.55

#81 Modified Apexo straight  
3.0 mm blade width.  
MDE.RE.56

Elevator #80 Straight gouge,  
3.0 mm blade width.  
MDE.RE.57

Heidbrink #40 Straight gouge,  
3.0 mm blade width.  
MDE.RE.58



Gentle back curve to improve adaptation



Heidbrink #41 Straight gouge  
3.7 mm blade width.  
MDE.RE.59

Heidbrink #1 Straight gouge  
4.5 mm blade width.  
MDE.RE.60

Stout #11A Gouge with  
tapered blade.  
MDE.RE.61

Fig. 11  
Stout  
MDE.RE.61X

Seldin #34,  
MDE.RE.62

Seldin #34S,  
MDE.RE.63



# Root Elevators | Root-Aufzüge

## Ascenseurs de racine | Elevadores de raíz



Potts #6 With hand form handle - left. MDE.RE.64	Potts #7 With hand form handle - right. MDE.RE.65	Crane Pick #8 Four-sided blade used for segmenting teeth MDE.RE.66	Hourigan #2 Left. MDE.RE.67	Hourigan #3 Right. MDE.RE.68
--	---	---	-----------------------------------	------------------------------------



Seldin #1L Left MDE.RE.69	Seldin #1R Right MDE.RE.70	Seldin #4L Left MDE.RE.71	Seldin #4R Right MDE.RE.72	Fig. 11M MDE.RE.731	Fig. 12M MDE.RE.732	Fig. 13M MDE.RE.733	Fig. 31F MDE.RE.741	Fig. 32F MDE.RE.742	Fig. 51 MDE.RE.751	Fig. 52 MDE.RE.752
---------------------------------	----------------------------------	---------------------------------	----------------------------------	------------------------	------------------------	------------------------	------------------------	------------------------	-----------------------	-----------------------



Cryer #25 Large flag - left. MDE.RE.76	Cryer #26 Large flag - right. MDE.RE.77	Cryer #44 Medium flag - left. MDE.RE.78	Cryer #45 Medium flag - right. MDE.RE.79	Cryer #64 Small flag - left. MDE.RE.80	Cryer #65 Small flag - right. MDE.RE.81
--	---	---	--	--	---




Miller Apexo #4 Left MDE.RE.82	Miller Apexo #5 Right MDE.RE.83	Modified Woodward #190 Left MDE.RE.84	Modified Woodward #191 Right MDE.RE.85	Woodward #15 Left MDE.RE.86	Woodward #16 Right MDE.RE.87
--------------------------------------	---------------------------------------	---	--	-----------------------------------	------------------------------------

# Root Elevators Bein | Root-Aufzüge Bein



Ascenseurs de racine Bein | Elevadores de raíz Bein

The 10° back curve improves adaptation.



Left & Right shank and compound blade angle allows easy access to 3rd molars.

Miller #71 Left MDE.RE.88  
Miller #72 Right MDE.RE.89



Miller #73 Left MDE.RE.90



Miller #74 Right MDE.RE.91



Coupland #1, 3.0 mm width and a square tip. MDE.RE.92



Coupland #2, 3.5 mm width and a square tip. MDE.RE.93



Coupland #3, 4.0 mm width and a square tip. MDE.RE.94

"English" pattern 2.5mm



Warwick James Straight MDE.RE.95  
Warwick James Right MDE.RE.96  
Warwick James Left MDE.RE.97



4.5mm 3.2mm 3.0mm  
1 HYLIN MDE.RE.98  
2 HYLIN MDE.RE.99A MDE.RE.99B



4.0 mm  
COLEMAN / LONDON HOSPITAL MDE.RE.111 MDE.RE.112 MDE.RE.113

7.5 mm

10 mm

13 mm

8 mm

13 mm



BARRY

1 MDE.RE.114 2 MDE.RE.115 3 MDE.RE.116 4 MDE.RE.117 17 MDE.RE.118 18 MDE.RE.119



WINTER

11L MDE.RE.120 11R MDE.RE.121 14L MDE.RE.122 14R MDE.RE.123



7 mm

1/1

LECLUSE MDE.RE.124



14mm Fig. 1 MDE.RE.125



14mm Fig. 2 MDE.RE.126



Fig. 1 MDE.RE.127



Fig. 2 MDE.RE.128



Fig. 4 MDE.RE.129




Fig. 5 MDE.RE.130





Root Elevators | Root-Aufzüge  
Ascenseurs de racine | Elevadores de raíz



*Original*

Root elevators Set  
Elévateurs  
Elevadores para raíces



Set of 8 Pieces  
MDE.RE.152X





*Original Bein*

4,5mm 3,5mm 2,5mm



BEIN 1 MDE.RE.145  
BEIN 2 MDE.RE.146  
BEIN 3

4,5mm 3,5mm



1# with angle MDE.RE.148  
2# MDE.RE.149

3,5mm 4,5mm



3# with acute angle MDE.RE.150  
4# MDE.RE.151



Popular

Straight Flared-Tip Elevators



ELEVATOR #52,  
2.75MM STRAIGHT,  
MDE.RE.160



ELEVATOR #77S,  
3.75MM STRAIGHT,  
MDE.RE.161



SERRATED  
Posterior Flared-Tip Elevators



ELEVATOR #77,  
3.75MM DISTAL ANGLED,  
MDE.RE.162



ELEVATOR #77R,  
3.75MM MESIAL ANGLED,  
MDE.RE.163

2,5mm 3,5mm 4,5mm 5,5mm 2,5mm 3,5mm 4,5mm 5,5mm



6 MDE.RE.164  
7 MDE.RE.165  
8 MDE.RE.166  
9 MDE.RE.167  
11 MDE.RE.168  
12 MDE.RE.169  
13 MDE.RE.170  
14 MDE.RE.171

Turning chisels



POTT



R MDE.RE.172  
L MDE.RE.173



Modell WIEN

Bajonett



3mm MDE.RE.174



4mm MDE.RE.175

Straight




3mm MDE.RE.176



4mm MDE.RE.177

Popular



Handle for interchangeable blades



**MDE.100.150.153**

Syndesmotom tips with handle, in case Contents: 20 tips and handle as depicted

**MDE.100.AN**

**Blue Line Extraction Elevator Kit**

5-piece extraction kit has Blue Titanium fortified handles with diamond coated tips.

: Straight, In, Out, Right, Left

*Diamond Coated*

**CHOMPRET**



1    2    3  
 MDE.RE.181    MDE.RE.182    MDE.RE.183



# Beck Modified Elevators



## Beck Modified Elevators



### TRI-BEVEL BACK SURFACE AND WAVE EDGE

Allows engagement and lift for the first 45° of elevation, at which point the edge will lose purchase from the tooth. Each bevel on the back of the Beck Elevators acts as its own fulcrum point, allowing the Beck elevators continuous engagement past the first 45° of rotation. This special sharpening, coupled with a thin tip for access into tight interproximal areas, make the Beck Elevators extremely effective instruments for elevation, but they must only be used by rotating, not prying them.

- MDE.100.192S**  
2.8MM, STRAIGHT
- MDE.100.191**  
3.5MM STRAIGHT
- MDE.100.192**  
4MM STRAIGHT
- MDE.100.191R**  
3.5MM, RIGHT
- MDE.100.191L**  
3.5MM, LEFT

### DOUBLE SIDED

Bernard For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.

Twisters for easily engage and cut periodontal ligaments. Featuring an unique twist to the blade that follows the contour of the root. The set of four contains two straight periostomes for anterior and two curved periostomes for the posterior

## Power Twist Elevators

### Set of 5

**MDE.99.AM**



- TW STR-L**  
ROOT TIP TWIST
- TW STR-R**  
ROOT TIP TWIST
- TW L**  
ROOT TIP TWIST
- TW R**  
ROOT TIP TWIST
- BERNARD 4RD**  
DOUBLE SIDED FACE  
GOLD PLATED TIP

### Set of 5 Bits and Handle **MDE.100.AM**



AM # 0



- AM # 1**  
Ø 1.6mm
  - AM # 2**  
Ø 2.8mm
  - AM # 3**  
Ø 3.4mm
  - AM # 4**  
Ø 2.8mm
  - AM # 5**  
Ø 2.8mm
- Working Length 15 mm

**ROCKMED Angle Modulation System**  
Minimally invasive widening of the alveolar ridge

# PEDODONTIC TOOTH FORCEP


## PEDIATRIC PATTERNS

### Child-Friendly Extractions

ROCKMED Rainbow colors to reduce child anxiety and create a more positive patient experience.



**ROCKMED**  
**Titanium Rainbow**



**MDE.DES.06**  
**Set 7 pcs Pedo Extraction Forceps Set**  
Contains: 2800/101-111-115-150-160-170-122  
Set children tooth forceps with spring

2800/Set 7 pcs PEDO  
Contains: 2800/101-111-115-150-160-170-122  
Set children tooth forceps with spring



2800/150 (mm 110)  
Lower incisors and canines

2800/170 (mm 110)  
Lower roots

2800/160 (mm 110)  
Lower molars and premolars



2800/101 (mm 110)  
Upper incisors



2800/111 (mm 110)  
Upper premolars



2800/115 (mm 120)  
Upper molars



2800/122 (mm 120)  
Upper roots

ROCKMED SQUARED BEAKS PATTERN



2400/133 BC Pedo (mm 130)  
Lower teeth, for children



2400/151 BC Pedo (mm 150)  
Lower teeth, for children



2400/101 BC Pedo (mm 150)  
Upper teeth, for children

# Extraction Forceps | Extraktion-Pinzette Daviers | Fórceps de extracción



## ENGLISH PATTERN *Child-Friendly Extractions*



2900/138

Upper teeth and roots, for children



2900/163

Upper centrals, for children



2900/37

Upper incisors, for children



2900/29 S

Upper roots, for children



2900/159

Upper premolars, for children



2900/157

Upper molars, for children



2900/158

Upper premolars and molars, for children



2900/39

Upper molars, for children



2900/39 A

Upper molars, for children



2900/39 R

Upper molars right, for children



2900/39 L

Upper molars left, for children



2900/38

Lower incisors and canines, for children



2900/162

Lower roots and centrals, for children



2900/123

Lower teeth and roots, for children



2900/13 S

Lower premolars, for children



2900/22 S

Lower molars, for children



2900/160

Lower molars, for children



2900/161

Lower molars, for children



2900/40

Lower molars, for children



# Extraction Forceps | Extraktion-Pinzette Daviers | Fórceps de extracción

Set Pedo tooth forceps, American pattern. 135mm  
3000/Set 9 pcs PEDO Contains:  
3000/65-69-16-17-18R-18L-203-150-151



AMERICAN PATTERN  
Amerikanische Form  
Forme américaine  
Forma americana  
Forma americana

Lower premolars, incisors, and roots



3000/151 Pedo 135mm  
(Cryer) universal  
3000/151S (165mm)



3000/203 Pedo 135mm  
Either side  
3000/203S (165mm)



3000/18 R Pedo (Harris, 135mm)  
Upper molars, right  
3000/18 RS (165mm)



3000/18 L Pedo (Harris, 135mm)  
Upper molars, left  
3000/18 LS (165mm)



3000/69 Pedo (Tomes, 135mm)  
Fragments, upper and lower  
3000/69S (165mm)



3000/16 Pedo (Harris, 135mm)  
Lower molars, cow-horn beak, either side  
3000/16S (165mm)



3000/17 Pedo (135mm)  
Lower molars, either side  
3000/17S (165mm)



3000/65 Pedo (135mm)  
Upper incisors and roots  
3000/65S (165mm)



3000/150 Pedo (Cryer, 135mm)  
Upper premolars, incisors, roots, universal  
3000/150S (165mm)

## Pedodontic tooth tweezers

Albrecht mm160

MDE.PTT.01



MDE.PTT.02



ROCKMED Rainbow colors to reduce child anxiety and create a more positive patient experience.



150SR\*  
Pedo



151SR\*  
Pedo

ROCKMED Rainbow

MDE.RB.150  
Upper Primary  
Teeth & Roots;  
Universal

ROCKMED Rainbow

MDE.RB.151  
Lower Primary  
Teeth & Roots;  
Universal

## Child-Friendly Extractions

Beak Modified Extraction Forceps, set of 7  
MDE.2020.7C

### ROCKMED Titanium Rainbow



MDE.RB.52  
# 52



MDE.RB.29  
# 29



MDE.RB.07  
# 7



MDE.RB.33  
# 33



MDE.RB.13  
# 13



MDE.RB.22  
# 22



MDE.RB.51  
# 51

# Signature Series

Refreshingly Innovative

Extraction Forceps | Extraktion-Pinzette | Daviers | Fórceps de extracción

## 2000 ROCKMED Signature Series

- anatomically correct
- prevent root fracture
- parallel gripping beak
- super heavy edge serrations
- available in 10 forceps designs
- contact both the crown and root structure
- long axis serrations within the beaks
- sharpened beaks to improve subgingival access
- used on intact, partially intact or broken down teeth



**ROCKMED CE**

Essential Extractions

Set of 10 CE

MDE.2020.10

2000 Model Series:

74, 02, 07, 21, 67, 79, 22, 17, 18, 13



ART# 45Z

Relief-Line

ART# 51Z

## Signature Series

### Noir 2000

Set of 10

MDE.2020.10N

Model Series:

74, 02, 07, 21, 67, 79, 22, 17, 18, 13

Finished in our new, non-reflective and incredibly wear resistant coating for a lifetime of function.



Premium Quality Dental Instruments  
Coidad premium dentales instrumentos de mano



LOWER ANTERIOR English Art#2074	UPPER ANTERIOR Art#2002	UPPER UNIVERSAL Art#2007	LOWER UNIVERSAL Art#2021	UPPER MOLAR UNIVERSAL Art#2067
---------------------------------------	----------------------------	-----------------------------	-----------------------------	--------------------------------------



LOWER ANTERIOR English Art#2013	UPPER RIGHT MOLAR Art#2018	UPPER LEFT MOLAR Art#2017	LOWER MOLAR English Art#2022	LOWER MOLAR Art#2079
---------------------------------------	----------------------------------	---------------------------------	------------------------------------	-------------------------





# Extraction Forceps | Extraktion-Pinzette Daviers | Fórceps de extracción

CERTIFIED COMPANY ISO

120 Series  
Standard Grip



140 Series  
Anatomical Grip



130 Series  
Diamond Coated



150 Series  
Feather-Light  
Grip



Diamond Coated Beak

Prefix 120, 130, 140, or 150 For Customized Grips



2100/1

Upper centrals and  
canines, wide



2100/2

Upper laterals and  
canines



2100/3

Upper crowded incisors  
and canines



2100/107

Upper canines



2100/180

Special for upper centrals  
and canines deep gripping



2100/181

Special for upper  
premolars deep gripping



2100/7

Upper premolars



2100/136

Upper centrals  
premolars and roots



2100/4

Lower incisors and canines



2100/75

Lower incisors  
and premolars



2100/75 A

Lower incisors  
and premolars



2100/137

Lower centrals and roots



2100/122

Lower centrals and roots



2100/182

Lower centrals canines and  
premolars deep gripping



2100/13

Lower premolars either side




2100/8

Lower premolars either side

# Extraction Forceps | Extraktion-Pinzette

## Daviers | Fórceps de extracción





# Extraction Forceps | Extraktion-Pinzette Daviers | Fórceps de extracción



2100/23 Lower molars, right    2100/24 Lower molars, left    2100/20 Lower third molars    2100/79 Lower third molars and lower wisdoms either side    2100/79 A Lower third molars    2100/56 For separating lower molars

## Cow horn beak



2100/86 Lower molars decayed or broken down crowns    2100/86 C Lower molars decayed or broken down crowns    2100/91 (Rowney) For separating lower molars    2100/99 (Wood) Lower molars    2100/87 Lower molars decayed or broken down crowns

## Upper roots



2100/51 Upper roots    2100/51 A Upper roots    2100/51 L Upper roots    2100/139 Upper roots    2100/52 Upper roots



2100/189 Fragments upper roots    2100/190 Fragments upper roots    2100/41 Upper roots    2100/49 Upper roots    2100/30 Upper roots

Extraction Forceps | Extraktion-Pinzette  
Daviers | Fórceps de extracción



2100/29

Upper roots and incisor



2100/113

Upper roots



2100/147

Upper roots



2100/101

Upper roots



2100/101 A

Upper roots



2100/108

Upper roots



2100/110 (Thorne)

Upper roots



2100/111

Upper roots



2100/76

Upper roots



2100/76 N

Upper roots



2100/76 S

Upper roots

Lower roots



2100/33

Lower roots



2100/33 A

Lower roots



2100/33 L

Lower roots



2100/74

Lower roots and incisors



2100/74 D

Lower roots



2100/74 M

Lower roots and incisors



2100/74 N

Lower roots and incisors



2100/114

Lower roots and centrals



2100/143

Lower roots and centrals

# Extraction Forceps | Extraktion-Pinzette Daviers | Fórceps de extracción



2100/187

Lower roots fragments



2100/46

Lower roots



2100/68

Lower roots



2100/185

Lower roots



2100/144 (Balding)

Lower roots and centrals



2500/350 (Routurier, mm 225)

Lower molars and wisdoms left



2500/351 (Routurier, mm 225)

Lower molars and wisdoms right



2500/130 (Ladmore, mm 190)

Upper wisdoms




2400/70

Foot fragments, horizontal embedded eights molars

## Extended "Deep Reach" Serrated Tips For Easy Atraumatic Extractions Tips Marked At 1mm-8mm

Deep Grip Extended Forceps  
Set of 4 Pieces: **MDE.DES.DG1**



Deep Gripping  
Upper Atraumatic Extraction Forcep  
**DGEF#A**



Deep Gripping  
Upper Atraumatic Extraction Forcep  
**DGEF#B**



Deep Gripping  
Lower Atraumatic Extraction Forcep  
**DGEF#C**



Deep Gripping  
Lower Atraumatic Extraction Forcep  
**DGEF#D**



UPPER UNIVERSAL  
**DGEF#E**



LOWER UNIVERSAL  
**DGEF#F**



UPPER ANTERIOR  
**DGEF#G**



UPPER 3RD MOLAR  
**DGEF#H**



LOWER MOLAR  
**DGEF#I**



LOWER ANTERIOR  
**DGEF#J**

Deep Grip  
ATRAUMATIC  
FORCEPS

Set of 6 Pieces:  
**MDE.DES.DG2**




MODIFY COW-HORN BEAK

# Extraction Forceps | Extraktion-Pinzette Daviers | Fórceps de extracción



for roots and fragments-TC



2500/221 (Witzel, mm 150) Universal for upper roots  
 2500/223 (Witzel, mm 150) Universal for lower roots  
 2175fig. 00 - LINDO LEVIEN Straight  
 2176fig. 01 - LINDO LEVIEN Curved  
 2177fig. 01 - WITZEL Lower roots  
 2178fig. 02 - WITZEL Upper roots



4650/24 Upper molars, either side  
 4650/18 R (Harris) Upper molars, right  
 4650/18 L (Harris) Upper molars, left  
 4650/88 L (Nevius) Upper molars, left  
 4650/88 R (Nevius) Upper molars, right



4650/10 S Upper molars either side  
 4650/210 Upper wisdoms, either side  
 4650/53 R Upper molars right  
 4650/53 L Upper molars left  
 4650/17 Lower molars, either side



4650/23 Lower molars, cow-horn beak, either side  
 4650/15 Lower molars, either side  
 4650/5 (Physick) Lower wisdoms, either side  
 4650/222 Lower wisdoms, either side



# Extraction Forceps | Extraktion-Pinzette

## Daviers | Fórceps de extracción



4650/150 (Cryer)

Upper incisors, premolars, roots, universal

4650/151 (Cryer)

Lower premolars, incisors, roots, universal

4650/150 A (Cryer)

Upper incisors, premolars, roots

4650/151 A (Cryer)

Lower premolars, incisors, roots

4650/150 AS

Upper premolars, incisors and roots

4650/151 AS

Lower premolars, incisors and roots



4650/32 A (Parmly)

Upper molars, premolars and roots, narrow beaks

4650/32 (Parmly)

Upper molars and premolars, broad beaks

4650/99 C

Upper premolars, canines and incisors

4650/65

Upper incisors and roots

4650/1 A

Incisor and canines

4650/286

Upper roots, premolars and incisors, either side



4650/203

Lower incisors, premolars and roots, either side

4650/103

Lower incisors, premolars and roots, either side

4650/85 A (Thomas)

Lower canines, premolars and molars, narrow beaks

4650/69 (Tomes)

Fragments, upper and lower

4650/6

For separating lower molars


4650/99 A (Kells)

Upper incisors, canines and premolars

### Tooth forceps scissor handle

For Lower wisdom teeth

MDE.TFS.88  
(mm 175)MDE.TFS.98  
(mm 175)



**MDE.DES.11** Complete Set of Apical Instrumentation - Apical Retention Forceps, PDL Luxators & Periostomes Forceps # 2, 7, 21, 67, 79, 74



**MDE.DES.01** / Set 10 pcs MEDFLAIR ESSENTIAL SET  
Contains: 2500/2-18A-51-22-75-33, elevators and periosteal



**MDE.DES.02** / Set 8 pcs MEDFLAIR ROOTS SET  
Contains: #51-189-29-46-187-74



**MDE.DES.03** / Set 10 pcs  
Contains: #22-13-79-187-74-1-7-17-18-51




**MDE.DES.05** / Set 12 pcs  
Contains: #2-13-33-51-7-79-17-18-22-67-74-86



**MDE.DES.04** / Set 10 pcs  
Contains: #2-13-33-51-7-79-17-18-22-67



**MDE.DES.08** / Set 10 pcs  
Contains: #2-13-33-51-7-79-17-18-22-67



**MEDFLAIR**  
Titanium Rainbow



**MDE.DES.06**  
Set 7 pcs Pedo Extraction Forceps Set  
Contains: 2800/101-111-115-150-160-170-122  
Set children tooth forceps with spring

**MDE.DES.07** / Set 10 pcs  
Contains: #2-13-33-51-7-79-17-18-22-67





Special Beak



**MDE.DES.12** / Modified Beak Forceps, Set of 8  
74, 02, 07, 21, 67, 79 & Spade Root Forceps



**MDE.DES.13** / Set 6 pcs ROCKMED Premium SET  
Extraction Forceps, Diamond Coated Beaks  
74, 02, 07, 21, 67, 79





# Atraumatic Forceps



150 Series  
Feather-Light  
Grip



**F1XS**

1 Standard Forceps  
Upper Incisors

**FAF1XS**

1 Apical Forceps  
Upper Incisors

**F67XS**

67 Forceps  
Upper Molars

**F35XS**

35 Forceps  
Upper Canines and Premolars

**FAF35XS**

35 Apical Forceps  
Upper Roots, Serrated

**FX51XS**

51 Root Forceps  
Upper Roots, Serrated

**F36XS**

36 Forceps  
Lower Incisors  
and Canines



**FX45XS**

45 Root Forceps  
Lower Roots, Serrated



**FX79XS**

79 Forceps  
Lower Molars



**FAF36XS**

36 Apical Forceps  
Lower Premolars, Canines  
and Incisors



**F13XS**

13 Forceps  
Lower Premolars



**FX73XS**

73 Forceps  
Lower Molars



**F23XS**

23 Cowhorn Forceps  
Lower Molars



ATRAUMAIR FORCEPS set of 13

MDE.DES.15

# "Beak and Bumper" Technique Extraction Forceps

Refreshingly Innovative



**ROCKMED CE**

Essential Extractions

Set of 4

Extraction Forceps | Extraktion-Pinzette | Daviers | Fórceps de extracción



Upper Anterior Forceps



Upper Right Forceps



Upper Left Forceps



Lower Universal Forceps

## Set code # MDE.2030

### Standard Series Set

- MDE.2030.1 Upper Right (Teeth #2-5)
- MDE.2030.2 Upper Left (Teeth #12-15)
- MDE.2030.3 Upper Anterior (Teeth #6-11)
- MDE.2030.4 Lower Universal (Teeth #18-31)
- MDE.2030.5 48 Pcs Bumper Guards



## Set code # MDE.2031

### Pedodontic Series Set

30% smaller in size

- MDE.2031.1 Upper Right (Upper Molars, Upper Canines)
- MDE.2031.2 Upper Left (Upper Molars, Upper Canines)
- MDE.2031.3 Upper Anterior (Incisors)
- MDE.2031.4 Lower Universal (2nd Molar to 2nd Molar)
- MDE.2031.5 48 Pcs Bumper Guards

**"Beak and Bumper" Technique Extraction Forceps four unique ROCKMED Forceps that allow you to easily extract any tooth.**

## Set code # MDE.2032

### Molar Series Set

- MDE.2032.1 Upper Left / Lower Right (Buccal)
- MDE.2032.2 Upper Right / Lower Left (Buccal)
- MDE.2032.3 24 Bumper Guards



## Set code # MDE.2033

FULCRUM EXTRACTION FORCEPS  
( Left and Right, set of 2 )



*Higher Flexibility*

**Ti-N**  
Titanium Nitride Coated

*Non-Stick*

FLEXIBLY RESILIENT AND SERRATED PERIOTOMES

**PT1SR** Double-Ended Flexible Serrated Anterior Periostome



Serrated instrument Designed to facilitate the removal of teeth with minimal damage to the surrounding alveolar bone by relieving the periodontal ligament. Ideal for anterior teeth, this instrument is particularly important when considering the placement of a dental implant.

**PT3SR** Double-Ended Flexible Serrated Posterior Periostome



Serrated instrument Designed to facilitate the removal of teeth with minimal damage to the surrounding alveolar bone by relieving the periodontal ligament. Ideal for posterior teeth, this instrument is particularly important when considering the placement of a dental implant.

Flexible - Two different flexibilities to address every Periostome challenge, Fracture resistant, Coated with titanium nitride for greater hardness and longer service life

**MDP.FLX3**  
Set of 2

**Self-Retaining Retractor**  
MICRO RETRACTOR FOR EDGES

Features:

- folding arms up to 40 degrees
- auto-retentive
- 3x3 sharp
- length 4.5 cm
- made of surgical steel



MDNT.0355

MODIFIED Agricola Retractor, 1.7in, 3X3 Prongs  
micro-dissecting retractors

*Easy - Out Periostomes*

Set of 7 Pieces  
**MDE.RE.PT10**

Periostomes are to be used as periodontal ligament knives that have fine tapering blades that compress the alveolar, cuts the membrane and gently eases the tooth from its socket. They were designed to make extractions as trauma free as possible. They are not to be used as elevators.

Backward Curve, 5mm    Straight, 5mm    Offset, 3mm    Forward Curve, 3mm    Backward Curve, 3mm    Straight, 3mm    Straight 2mm



# Handle Designs ROCKMED®

Refreshingly Innovative

**SINOX**  
Stainless Steel

The Handle Versions

*Safe and easy to clean*

Standard

**A**



**B**



**C**



Bein

**D**



**E**



Finger Form

**F**



**G**



Carpal

**H**



**J**



The ROCKMED Logo, and the Name are Registered and Trade Mark of ROCKMED Industries



Scan our mobile website to your smartphone for quick reference



## Corporate Head

61-B Liaqat Ali Road, Industrial Estate,  
Sialkot, Pakistan

Phone # +92-523-555526

Fax # +92-523-556226

+ 92-333-8623458

+92 331 3177777

info@rockmed.co

www.rockmed.co



*Precision Craftsmanship*

CERTIFIED COMPANY  
ISO-9001 ISO-13485



Copyright © ROCKMED Industries  
Reprinting of this catalog even partially is not allowed.  
ROCKMED Management Information System.